

## Introduction

The StarGazer-2 (*Epiphyte3*) instrument is designed for low-volume, 384-well optical bottom microplates used for high throughput data collection.




- Label-free thermo-aggregation assay using differential static light scattering (DSLS) at 620 nm.

## Protein Applications of the StarGazer2:

Technical Specifications	
Sample volume	2.5 - 50 $\mu$ L
Temperature range	10 - 95 $^{\circ}$ C
Heating rate	0.1 - 5 $^{\circ}$ C/minute
Protein conc. range	0.05 - 1 mg/mL
DSLS laser wavelengths	620 nm
DSLS resolution	$\sim$ 8 kDa mean
Number of Samples	1 - 383

- Temperature gradient and isothermal-based experiments
- Optimal buffer screening
- Compound stabilization / destabilization screening and dose response characterization
- Nucleic acid and peptide binding characterization
- Evaluation of protein refolding conditions
- Chemical fingerprinting
- Formulation development of therapeutic monoclonal antibodies
- Comparison of stability of SNP proteins



**Structural & Biophysical Core Facility**

## **Greg Wasney**

**Manager, Structural & Biophysical Core Facility**  
Peter Gilgan Centre for Research & Learning  
The Hospital for Sick Children  
686 Bay Street, Rm. 21.9708  
Toronto, ON. M5G 0A4

Email: [greg.wasney@sickkids.ca](mailto:greg.wasney@sickkids.ca)

Office: 416.813.7209  
Office Internal Ext. 307209  
Lab: 416-813-7654 ext. 309442

<http://lab.research.sickkids.ca/sbc-facility/>

**SickKids**  
THE HOSPITAL FOR  
SICK CHILDREN

## The Hospital for Sick Children's Structural & Biophysical Core Facility



## StarGazer-2— Epiphyte3

Differential Static Light Scattering (DSLS) Protein Thermo-denaturation

## Application Note #1:

Optimal protein buffer screen

- Greater long/short-term stability
- Reach higher concentration, recover low solubility proteins
- Better binding / functional activity
- Promotion of protein crystallization

# StarGazer-2— Protein Thermodenaturation Assays

## Experiment Type 1: Temperature



### Gradient Denaturation

**Protein Sample:** 0.5 mg/mL, 10  $\mu$ L

**Temperature ramp rate:** 1  $^{\circ}$ C/min, 25  $^{\circ}$ C to 95  $^{\circ}$ C

**Buffers:** 100 mM MES pH 6.0; 100 mM HEPES pH 7.0, 8.0; 100 mM Glycine pH 9.0

**Salt:** 50, 250, 500 mM NaCl



Best: Low pH, High [NaCl]

The protein is stabilized by buffers with lower pH and higher NaCl concentrations.

## Experiment Type 2: Isothermal Denaturation

**Protein Sample:** 0.5 mg/mL, 10  $\mu$ L


**Isothermal denaturation @ 37  $^{\circ}$ C for 200 min**

**Buffers:** 100 mM MES pH 6.0; 100 mM HEPES pH 7.0, 8.0; 100 mM Glycine pH 9.0

**Salt:** 50, 250, 500 mM NaCl

**Sypro Orange:** 5X

Performed initial temperature gradient experiment using DSF (differential scanning fluorometry) to determine temperature of denaturation using a midrange buffer of 100 mM HEPES, pH 7.5, 250 mM NaCl.




Based on the protein thermo-denaturation gradient curve above, an isothermal denaturation temperature of 37  $^{\circ}$ C should be used (2  $^{\circ}$ C below the onset of thermost-denaturation at 39  $^{\circ}$ C).


Therefore, the isothermal temperature was held at 37  $^{\circ}$ C for 200 minutes.

The protein is stabilized by buffers with higher NaCl concentrations while no significant changes are observed with varying pH.

### Fast Denaturation (50mM NaCl)



### Slower Denaturation (250mM NaCl)



### Little Denaturation (500mM NaCl)

