

Introduction

The StarGazer-2 (*Epiphyte3*) instrument is designed for low-volume, 384-well optical bottom microplates used for high throughput data collection.



- Label-free thermo-aggregation assay using differential static light scattering (DSLS) at 620 nm.

Protein Applications of the StarGazer-2:

Technical Specifications	
Sample volume	2.5 - 50 μ L
Temperature range	10 - 95 $^{\circ}$ C
Heating rate	0.1 - 5 $^{\circ}$ C/minute
Protein conc. range	0.05 - 1 mg/mL
DSLS laser wavelengths	620 nm
DSLS resolution	\sim 8 kDa mean
Number of Samples	1 - 383

- Temperature gradient and isothermal-based experiments
- Optimal buffer screening
- Compound stabilization / destabilization screening and dose response characterization
- Nucleic acid and peptide binding characterization
- Evaluation of protein refolding conditions
- Chemical fingerprinting
- Formulation development of therapeutic monoclonal antibodies
- Comparison of stability of SNP proteins



Structural & Biophysical Core Facility

Greg Wasney

Manager, Structural & Biophysical Core Facility
Peter Gilgan Centre for Research & Learning
The Hospital for Sick Children
686 Bay Street, Rm. 21.9708
Toronto, ON. M5G 0A4

Email: greg.wasney@sickkids.ca

Office: 416.813.7209
Office Internal Ext. 307209
Lab: 416-813-7654 ext. 309442

<http://lab.research.sickkids.ca/sbc-facility/>

SickKids
THE HOSPITAL FOR
SICK CHILDREN

The Hospital for Sick Children's Structural & Biophysical Core Facility



StarGazer-2 — Epiphyte3

Differential Static Light Scattering (DSLS) Protein Thermo-denaturation

Application Note #3:

DNA and peptide stabilization

- Protein-DNA interaction mapping
- Protein-Peptide interaction mapping
- Protein-Polymer interaction mapping

StarGazer-2— Protein Thermo-denaturation Assays

Protein-DNA interaction mapping

Protein Sample: 25 μ M, 10 μ L

Temperature ramp rate: 1 $^{\circ}$ C/min, 25 $^{\circ}$ C to 95 $^{\circ}$ C

Buffer: 100 mM HEPES pH 8.0, 250 mM NaCl

DNA concentration: 50 μ M

Performed temperature gradient experiments to determine if protein stability is affected by the presence of specifically sized DNA.



Addition of a full length DNA fragment significantly increased the protein stability due to protein-DNA interaction.

Addition of a largely truncated fragment (Tr.) only slightly increases protein stability. However, as the truncated fragment was elongated by two (Tr. +2) or four (Tr. +4) base pairs, its stabilizing effect on the protein increases.

The truncated fragment elongated by 4 base pairs had a similar stabilizing effect as the full length fragment. Thus, the Protein-DNA interaction can be specific to a region approximately within the Tr. +4 fragment.

This protein thermo-denaturation based interaction mapping strategy can be similarly applied to peptide and polymer binding.

DNA Fragment	Stability
No DNA	-
FL	+++
Tr.	+
Tr. +2	++
Tr. +4	+++



Protein stability is increased in the presence of specifically sized DNA.