40 ul of pGIPZ non-silencing supernatant ~6.5x10e6 TU/mL added to 24 well plate with 225ul media with MDCK-75000 cells/well: 48 hour incubation (GFP image)

40 ul of pGIPZ non-silencing supernatant ~6.5x10e6 TU/mL added to 24 well plate with 225ul media with MDCK-75000 cells/well: 48 hour incubation (merged image)

60 ul of pGIPZ non-silencing supernatant ~6.5x10e6 TU/mL added to 24 well plate with 225ul media with MDCK-75000 cells/well: 48 hour incubation (GFP image)

60 ul of pGIPZ non-silencing supernatant ~6.5x10e6 TU/mL added to 24 well plate with 225ul media with MDCK-75000 cells/well: 48 hour incubation (merged image)